


1. YIWEI JIA ON USE-CASE 2: EYE-TRACKER



Use case 2: Eye-Tracker



Questions:


- Channel: Human eye-movements via infra-red camera
- Environment: Dark room, possible head movements
 - Noise: Blinks, Eye-movements, Sampling failure (device error)
 - How can I filter this noise: Erase data points affected by artefacts
- Dimensionality of the signal: 1D and 2D (why?)
- What is the sampling rate? Infra-red camera specifics
- Sample size: 4

Which process am I investigating?

- Gaze stability (fixation ability of the subjects), Deviation from center






Use case 2: Eye-Tracker




Context of eye tracking data processing:

- **Motion artifact** : slow exposure + moving subjects = **blur**.




Use case 2: Eye-Tracker




Context of eye tracking data processing: Motion artifact

Long acquisition time of MRI +

Head movement




Rapid eye movement




Use case 2: Eye-Tracker

Context of eye tracking data processing:


- **Motion artifact :**
Long acquisition time of MRI + **Rapid** eye movement
- **Method: Track the eye movement in real-time**
 - Records: gaze position
 - Types of eye movement: *Blinks, Saccades, Tremor*
 - Goal: *eliminate MRI data affected by motion*



A dot at the center of screen



Gaze positions (one MRI session)




Use case 2: Eye-Tracker

Context of eye tracking data processing:


- **Motion artifact :**
Long acquisition time of MRI + **Rapid** eye movement
- **Method: Track the eye movement in real-time**
 - Records: gaze position
 - Types of eye movement: *Blinks, Saccades, Tremor*
 - Goal: *eliminate MRI data affected by motion*



A dot at the center of screen



Gaze positions (one MRI session)




Use case 2: Eye-Tracker

Context of eye tracking data processing:

- **Motion artifact :**
Long acquisition time of MRI + **Rapid** eye movement
- **Method: Track the eye movement in real-time**
 - Records: gaze position
 - Types of eye movement: **Blinks, Saccades, Tremor**
 - Goal: **eliminate data affected by motion**
- **What we do with ET data**
 - Filter Noise: device error (eg. extremely large values)
 - Eliminate data points affected by motions.
 - Investigate Eye Movement Features
 - Understand the characteristics, quality control

Use case 2: Eye-Tracker




Part 1: Know your data at the first glance (parsing metadata/ visualization)

- Load the data and parse the metadata
- Plot each dimension of gaze trajectories over time
- Overlay different event types (saccades, blinks, fixations) using different colors or markers.

What is the correlation between the events and gaze trajectories?


- Do **fixations** correspond to relatively **flat** or stable gaze lines?
- Do **saccades** appear as **sharp jumps**?
- Are there any periods where the **gaze suddenly drops out**, possibly indicating a **blink**?

Use case 2: Eye-Tracker





Take Subject 1 as example:

Check the gaze trajectory in two dimensions along time axis




Use case 2: Eye-Tracker

Check the **pupil area** in two dimensions along time axis




Zoom in

Zoom In - Pupil Area




discontinuous
Profiles ?

Use case 2: Eye-Tracker





Overlap the **blink event** with pupil area





Use case 2: Eye-Tracker

Check the x coordinate of gaze trajectory with **saccade**




Use case 2: Eye-Tracker

Visualize the gaze point in 2D



Points deviated from the center of the screen

Use case 2: Eye-Tracker






Part 2: Feature analysis


- **Distribution** of gaze points along each direction.
 - Determine the **center** of gaze points and the **deviation**.
 - Calculate **statistics** (mean, median, standard deviation) for each subject
 - Saccade duration, blink duration, and fixation duration.
- Solid set of features per subject
- Quality analysis
 - Do the features make sense compared to expected behavior?
 - Do some subjects need to be excluded from further analysis due to poor data quality?

Use case 2: Eye-Tracker

Filter out the gaze points affected by **blinks** and **saccades**
Preserve data during the **fixation**



Use case 2: Eye-Tracker



Calculate the statistics of the gaze points along each dimension

Statistic	x-Dimension	y-Dimension
Mean	397.07	305.27
Median	395.5	305.30
Standard Deviation	13.96	9.86
Minimum	353.40	225.80
Maximum	609.80	460.70
25th Percentile (Q1)	388.30	301.00
75th Percentile (Q3)	405.30	309.00


Recall: Cours4

Descriptive statistics of a time series

Signal characterization:


1. **Max, Min, Mean** - Central value of the signal.
2. **Median**: Middle value when sorted.
3. **Variance & Standard Deviation**: Measure of signal spread (and fluctuation).
4. **Distribution** analysis: Shows patterns like normal or skewed distribution.

Use case 2: Eye-Tracker




Investigate the density of gaze points along each dimension (before filtering vs. after filtering)

X Coordinate: Raw vs. Filtered




Y Coordinate: Raw vs. Filtered



How has the distribution changed after filtering?

Use case 2: Eye-Tracker




Check Quartiles: Values that divide a data set into four equal parts after it has been sorted.

- Q1 (25th percentile):** The value that separates the lowest 25% of the data. 25% of the data points are less than Q1.
- Q2 (median):** The middle value of the data set, which divides the data into two equal halves.
- Q3 (75th percentile):** The value that separates the lowest 75% of the data. 75% of the data points are less than Q3.

IQR: the difference between the 75th percentile (Q3) and the 25th percentile (Q1)

LIBRE




Upper boundary = Q3


Lower boundary = Q1

Small IQR -> The majority of gaze points are close to the median



LIBRE




Use case 2: Eye-Tracker





Gaze points located between Q1 and Q3




Use case 2: Eye-Tracker





Adjust different upper and lower boundaries according to specific needs



Use case 2: Eye-Tracker




Corresponding filtering results



Use case 2: Eye-Tracker

Investigate the eye movement features across 4 subjects



Blink (ms)

Subject	Mean Duration	Median Duration	Standard Deviation	Min Duration	Max Duration	Total Duration	Number of Times
1	307.87	298.0	102.14	136	612	9544	31
2	193.79	191.0	53.48	84	366	8333	43
3	279.44	257.0	95.56	82	643	17605	63
4	570.57	222.5	963.50	2	7006	94714	166

Saccade (ms)


Subject	Mean Duration	Median Duration	Standard Deviation	Min Duration	Max Duration	Total Duration	Number of Times
1	30.50	17.0	63.47	6	612	22203	728
2	36.31	22.0	50.18	6	366	19206	529
3	51.54	26.0	79.13	6	643	36796	714
4	151.31	63.0	390.59	2	7006	202600	1339

Subject 004
looks different?

blinks and saccade
→ Long total duration
→ Occurred many times

Use case 2: Eye-Tracker

Investigate the eye movement features across 4 subjects



Fixation

Subject	Mean Duration	Median Duration	Standard Deviation	Min Duration	Max Duration	Total Duration	Number of Times
1	485.33	247.0	624.00	22	4157	353804	729
2	672.02	491.5	594.36	17	5003	356169	530
3	417.98	246.0	466.94	24	3348	298856	715
4	128.82	92.0	117.28	9	1012	172486	1339

Subject 004
looks different?


Fixation

→ Short duration

Reflect the data quality




✓ Think about quality control

Use case 2: Eye-Tracker



Visualize the statistics of event durations across four subjects.

Violin plots are handy tools



Think about quality control of early data:

- Poor eye tracker calibration?
- Subject was tired?
- ...

Is good for



- ✓ Retrospectively update the SOP of your experiment.
- ✓ Predict results in future steps.
- ✓ Stop the processing loop and save more time

Hand-on session: Eye-Tracker

https://github.com/MattechLab/SigInf_2025

<https://colab.research.google.com/>

